

INDIAN INSTITUTE
OF

TECHNOLOGY

PHYSICS DEPT.

NO.

NAME A-SRIYA

EXPT. No. 8

CLASS

BATCH

LABORATORY

DATE 26-10-23

FRESNEL'S BI-PRISM

Aim: To obtain the wavelength of Na light using fresnel's biprism

Apparatus : (i) Optical rail with uprights

(ii) Sodium lamp


(iii) Bi-prism

(iv) Convex lens

(v) Slit & micrometer-eyepiece

Theory :

- 1) The fresnel's bi-prism has one angle slightly less than two right angles and two equal, small base angles act like 2 very thin prism placed base to base
- 2) When the rays from slit S illuminated by mono-chromatic light like Na light are incident on bi-prism, The emergent rays from 2 halves seems (appears) to diverge from 2 virtual sources S_1 & S_2 .



- 3) If the screen is placed in plane 1 to slits, emergent rays are in to overlap to form dark and bright fringes
- 4) If d (dis b/w S_1 & S_2); then relation is

$$\boxed{x = \frac{dD}{d}} \quad \text{--- (1)}$$



$$\boxed{d = \frac{d}{D} x} \quad \text{--- (2)}$$

D - dis b/w slit & screen
 d - wavelength of light
 x - fringe width

- 5) To determine d , convex lens is placed between slit and screen with focal length less than one fourth of the distance b/w slit & screen
- 6) lens is adjusted so that to obtain real images of the virtual sources S_1 & S_2 on the screen with d_1 & d_2 are their separations in 2 positions \Rightarrow

$$\boxed{d = \sqrt{d_1 d_2}} \quad \text{--- (3)}$$

Removal of lateral shift



INDIAN INSTITUTE
OF

TECHNOLOGY

PHYSICS DEPT.

NO.

NAME

EXPT. No.

DATE

CLASS

BATCH

LABORATORY

Formulae used:

$$1) d_{avg} = \frac{d_1 + d_2}{2}$$

$$\Delta d = \left| \frac{d_1 - d_2}{2} \right|$$

2) Fractional error in x :

$$\frac{\Delta x}{x} = \left| \frac{x' - x''}{2x} \right|$$

; x' & x'' are the slopes of x_n vs n graph

3) As $d = \frac{x d}{D}$

$$\frac{\Delta d}{d} = \frac{\Delta x}{x} + \frac{\Delta d}{d} \quad \left(\left(\frac{\Delta D}{D} \right) \rightarrow \text{negligible} \right)$$

Precautions:

1) Dis. b/w slit and biprism ^{or slit width} should be kept constant

2) Move the eye-piece only backward for adjusting lateral shift, move bi-prism laterally

3) During measurement of d_1 & d_2 ; Do not keep eye-piece at large distance from slit ($\gg 4f$) that results in great error

4) While measuring fringe width, align cross wire at bright fringe

5) While measuring fringe width, the cross-wire must be moved only in 1 dir. to avoid backlash

Observations

* Distance between slit and eyepiece $D = 1459\text{mm} - 459\text{mm}$

Measurement of fringe width. $\frac{1000}{1000} \text{ mm} = 100 \text{ cm}$

Least count

SNo of fringes	Main scale reading MSR = (mm)	Circular scale reading	Total reading $\times n$ (mm)
1)	14	30 \Rightarrow 0.30mm	14.30
2	13.5	46 \Rightarrow 0.46mm	13.96
3	13.5	16 \Rightarrow 0.16mm	13.66
4	13	35 \Rightarrow 0.35mm	13.35
5	13	5 \Rightarrow 0.05mm	13.05
6	12.5	33 \Rightarrow 0.33mm	12.83mm
7	12.5	4 \Rightarrow 0.04mm	12.54mm
8	12	27 \Rightarrow 0.27mm	12.27mm
9	11.5	47 \Rightarrow 0.47mm	11.97mm
10	11.5	20 \Rightarrow 0.20mm	11.70mm

Measurement of dis b/w virtual sources d :

$$d_1 = \frac{1459\text{mm}}{1000} \cdot 1.80\text{mm} = \frac{18.65 - 10.80}{1000} = 7.85\text{mm}$$

$$d_2 = \frac{750\text{mm}}{1000} \cdot 2.49\text{mm} = \frac{18.35 - 17.93}{1000} = 0.42\text{mm}$$

$$d = \sqrt{d_1 d_2} = \sqrt{7.85 \cdot 0.42} = 1.82\text{mm}$$

for calculation of error; measuring d at diff eyepiece posn

$$d_1' = \frac{1164\text{mm}}{1000} \cdot 20.03 - 13.71 = 6.32\text{mm}$$

$$d_2' = \frac{689\text{mm}}{1000} \cdot 17.73 - 16.79 = 0.94\text{mm}$$

$$d' = \sqrt{d_1' d_2'} = \sqrt{6.32 \cdot 0.94} = 2.43\text{mm}$$

INDIAN INSTITUTE
OF

TECHNOLOGY

PHYSICS DEPT.

NO.

NAME

EXPT. No.

DATE

CLASS

BATCH

LABORATORY

Calculations:

1) Fringe width x (from graph)

x' & x'' are max & min slopes of graph of the lines

$$x = \frac{x' + x''}{2} = \frac{0.3125 + 0.2800}{2} = 0.2962 \text{ mm}$$

$$\Delta x = \left| \frac{x'' - x'}{2} \right| = \frac{0.3125 - 0.2800}{2} = 0.0162 \text{ mm}$$

2) Virtual source distance d :

$$d = \frac{d + d'}{2} = \frac{4.253 \text{ mm} + 2.125 \text{ mm}}{2} = 2.125 \text{ mm}$$

$$\Delta d = \left| \frac{d' - d}{2} \right| = \frac{16.99 \text{ mm} - 2.125 \text{ mm}}{2} = 0.305 \text{ mm}$$

3) wavelength λ : ; $\lambda = \frac{d}{D}(x)$

$$= \left(\frac{2.125}{1000} \right) (0.2962) \times 10^3 \text{ m} = 6.29425 \times 10^{-7} \text{ m} = 629.425 \text{ nm}$$

$$\frac{\Delta \lambda}{\lambda} = \frac{\Delta x}{x} + \frac{\Delta d}{d}$$

$$= \left(\frac{0.0162}{0.2962} + \frac{0.305}{2.125} \right) = 0.05469 + 0.14353 = 0.198$$

$$\Delta \lambda = \left(\frac{\Delta \lambda}{\lambda} \right) (\lambda)$$

$$= (0.198)(629.425) \text{ nm}$$

$$= 124.76 \text{ nm}$$

RESULT:

Using the Fresnel's biprism experiment, the wavelength of Na light comes out to be

$$\lambda = (\lambda \pm \Delta \lambda) = (629.425 \pm 124.76) \text{ nm}$$

Graph of x_n (vs) n .

Scale:

x-axis, 1 unit = 1

y-axis, 1 unit = 0.2 mm

$$|\text{max slope } x'| = \left| \frac{14.3 - 13.0}{5 - 1} \right| = 0.3125 \text{ mm}$$

$$|\text{min slope}| = |x''| = \left| \frac{13.66 - 11.7}{10 - 3} \right| = 0.280 \text{ mm}$$

(10)

26

